

## Product datasheet for TA130001

## Goat Anti-Mouse IgG Secondary Antibody

## **Product data:**

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Secondary Antibodies
Product Name:	Goat Anti-Mouse IgG Secondary Antibody
Applications:	ELISA, ICC/IF, IHC, WB
Recommended Dilution:	WB, IHC, IF/ICC, ELISA
Reactivity:	Mouse IgG
Host:	Goat
Immunogen:	Mouse IgG, whole molecule
lsotype:	lgG
Formulation:	Lyophilized with 1% sodium acetate
Purification:	This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all goat serum proteins, except the specific antibody for mouse lgG.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US