

Product datasheet for **TA150102**

Clone OTI4G1, Anti-GST mouse monoclonal antibody


Product data:

Product Type:	Tag Antibodies
Product Name:	Clone OTI4G1, Anti-GST mouse monoclonal antibody
Clone Name:	OTI4G1 (formerly 4G1)
Applications:	WB
Recommended Dilution:	WB: 1:5000
Target:	GST tag
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	Recombinant GST produced in E.coil
Isotype:	IgG1
Formulation:	0.13mg/ml in PBS (pH7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide
Concentration:	0.13mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	Stable for 1 year at -20°C from date of shipment.
Synonyms:	ANTI-GST; ANTIBODY FOR GST; ANTI GST; GST AB; GST ANTIBODY



[View online »](#)

Product images:



Anti-GST immunoblot analysis of six different HEK293T over-expression lysates containing different GST tagged recombinant proteins, and one GST-DDK bacterial lysate produced in E.coli. at 1:5000 dilution.