

Mouse IL-1 β (IL-1 beta, IL-1F2) Protein, Recombinant

产品编号	产品名称	包装规格
NBS0256-5 μ g	Mouse IL-1 β (IL-1 beta, IL-1F2) Protein, Recombinant	5 μ g
NBS0256-10 μ g	Mouse IL-1 β (IL-1 beta, IL-1F2) Protein, Recombinant	10 μ g
NBS0256-20 μ g	Mouse IL-1 β (IL-1 beta, IL-1F2) Protein, Recombinant	20 μ g
NBS0256-50 μ g	Mouse IL-1 β (IL-1 beta, IL-1F2) Protein, Recombinant	50 μ g
NBS0256-100 μ g	Mouse IL-1 β (IL-1 beta, IL-1F2) Protein, Recombinant	100 μ g

产品简介:

IL-1 beta/IL-1F2 Protein, Mouse, Recombinant is expressed in E. coli expression system. The predicted molecular weight is 17.5 kDa and the accession number is P10749.

产品性质:

生物活性	Measured by its ability to induce Interferon gamma secretion by human natural killer lymphoma NK-92 cells. The ED50 for this effect is typically 6-30 ng/mL.
产品别名	白细胞介素, 白介素, interleukin 1, β , interleukin 1, beta, IL-1 β , IL-1beta, IL-1b, IL-1 β
种属	Mouse
蛋白编号	P10749
标签	No Tag
表达系统	E.coli
分子量	17.5 kDa (predicted); 19 kDa (reducing conditions)
纯度	> 97% as determined by SDS-PAGE
内毒素	< 1.0 EU per μ g by the LAL method.
缓冲液	Lyophilized from a solution filtered through a 0.22 μ m filter,

	containing 20 mM Tris, pH8.5. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.
外观	Lyophilized Powder
复溶方法	Reconstituted with sterile deionized water to 0.25 mg/mL.

保存条件:

冻干粉-20℃ ~ -80℃保存, 1 年有效; 复溶后, 无菌条件下, 2~8℃保存, 7 天有效。
复溶后, 无菌条件下, -85℃ ~ -65℃保存, 3 个月有效。避免反复冻融, 并将产品分装储存。

注意事项:

1. 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

本产品仅用于生命科学研究, 不得用于医学诊断及其他用途!